Prontosan® - Wound irrigation solution

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en Prontosan® - Wound irrigation solution

Intended Purpose: For rinsing, cleansing, moistening and decontamination of acute, chronic, and infected skin wounds and burns.

1. Introduction: Chronic skin wounds are often coated with slough, necrotic tissue and/or biofilm. These coatings are difficult to remove and lead to delayed wound healing. Prontosan® Wound Irrigation Solution is able to remove these barriers to wound healing through its cleansing activity. Acute wounds also require proper cleansing as they are generally contaminated with debris and microorganisms. These contaminants can interfere with the normal wound healing process and lead to complications such as infection. Due to the unique combination of ingredients (i.e. the antimicrobial substance, polyhexanide, and the surfactant, betaine), Prontosan® Wound Irrigation Solution is ideal for the prevention of biofilm formation. Prontosan® Wound Irrigation Solution should be used for wound cleansing prior to treatment with Prontosan® Wound Gel or Prontosan® Wound Gel X. The choice between Prontosan® Wound Gel and Gel X allows an optimal application in large surface areas as well as deep wounds. In general, both products are suitable in all conditions. However, due to the higher viscosity of Prontosan® Wound Gel X, we recommend to use it preferably for large surface area wounds including burns, whereas Prontosan® Wound Gel is more fluid and might therefore be easier to apply in deep or tunnelling wounds, wound cavities and difficult to access areas. The summary of safety and clinical performance (SSCP) is available in the European database on medical devices (Eudamed) and can be found using basic UDI-DI 4039239000009242S or the product name. <https://ec.europa.eu/tools/eudamed>. Or requested from the Manufacturer.

2. Product profile and areas for use: For rinsing, cleansing, moistening and decontamination of:

- acute non-infected and infected wounds; traumatic wounds (e.g. lacerations, crush injuries);
- chronic non-infected and infected wounds including pressure ulcers; vascular ulcers; diabetic ulcers.
- postoperative wounds.
- thermal and non-thermal burns degree: IIa, IIb and III.
- fistulas and abscesses.
- entry ports of urological catheters, PEG/PEJ tubes or drainage tubes.
- peristomal skin.

For intraoperative cleansing and irrigation of wounds. For moistening encrusted dressings and bandages prior to removal.

3. General use: To ensure optimal cleansing results with Prontosan® Wound Irrigation Solution, aseptic application is recommended whenever possible. Prontosan® Wound Irrigation Solution may be used to saturate compresses, gauzes or pads for cleansing as needed. Application should be conducted with every dressing change in order to achieve and maintain a visually clean wound. Prontosan® Wound Irrigation Solution may be warmed to body temperature prior to use. Wound dressings and bandages are frequently encrusted and adhere to the wound surface. Removal of these dressings can traumatize the wound bed, which in turn delays wound healing. In cases where dressings and bandages are difficult to release due to encrustation, soaking the dressing with Prontosan® Wound Irrigation Solution is recommended until they can be gently removed. Transurethral catheters: Thoroughly moisten a sterile gauze and cleanse the catheter entry site.

Peg- / Pej-tubes and suprapubic catheters: For cleansing of entry port, apply a sufficient amount of Prontosan® Wound Irrigation Solution. Encrustations can be moistened with Prontosan® Wound Irrigation Solution in order to be easily removed. Dry surrounding skin area with a sterile gauze and apply a sterile cover afterwards.

Stoma care: Thoroughly moisten a sufficient number of sterile compresses. Clean peristomal skin areas and dry afterwards. Prior to the use of Prontosan® Wound Irrigation Solution, rinse off the residues of soaps, ointments and oils or other wound cleansing substances thoroughly.

4. Intended users: The Prontosan® Wound Irrigation Solution is intended to be used by professional users or by instructed adults only.

The intended patient population are adults and children.

5. Wound cleansing: The entire wound, as well as the surrounding area, should be thoroughly cleansed with Prontosan® Wound Irrigation Solution. For non-severe and acute wounds rinsing is sufficient. For chronic wounds, soaking the affected area with Prontosan® Wound Irrigation Solution for at least 5 minutes with a saturated compress is recommended. The solution should be used undiluted.

6. Tissue tolerability and biocompatibility: Dermatologically tested and evaluated as non-irritating and well-tolerated; painless; no inhibition of granulation or epithelialization.

7. Side effects: In very rare cases, there may be a mild burning sensation after application of Prontosan® Wound Irrigation Solution, but this usually dissipates after a few minutes. Prontosan® Wound Irrigation Solution can cause allergic reactions such as itching (urticaria) and rashes (exanthema). In very rare cases (less than 1 out of 10,000), anaphylactic shock has been reported.

8. Contraindications: Prontosan® Wound Irrigation Solution should not be used:

- if the patient is known to be allergic or if it is suspected that the patient may be allergic to one of the ingredients of the product.
- on the central nervous system or the meninges.
- in the middle or inner ear.
- in the eyes. If Prontosan® Wound Irrigation Solution comes into contact with the eyes, flush the eyes with running water and seek medical advice.
- on hyaline cartilage and in aseptic joint surgery. If Prontosan® Wound Irrigation Solution does come into contact with aseptic cartilage, it should be immediately irrigated with Ringer's solution or normal saline.
- for peritoneal lavage/irrigation.

9. Restrictions of use: Pregnancy and lactation period: There is no evidence of mutagenic or embryo toxicity associated with the ingredients of this product. Due to the lack of relevant clinical trials and clinical experience with pregnant and breast-feeding women, Prontosan® Wound Irrigation Solution should not be used.

10. General safety instructions: Do not use for infusion, injection or ingestion. Do not use in combination with cleansing soaps (anionic surfactants), ointments, oils, etc.


These substances should be thoroughly removed from the wound before use as they can affect the cleansing performance of Prontosan® Wound Irrigation Solution. Do not use damaged bottles. Keep bottles away from direct sunlight. Keep out of reach of children. The use of Prontosan® Wound Irrigation Solution does not replace the need for systemic or other adequate infection treatment.

11. Summary / technical information: Prontosan® Wound Irrigation Solution is a preserved product and has a shelf life of 8 weeks after opening. The bottle should be closed immediately after use to prevent contamination. The top of the bottle should be protected from contamination during use. Bottles that have come into direct contact with the wound or have become contaminated in another way must be discarded. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Composition: Purified Water, Betaine Surfactant, 0.1 % Polyaminopropyl Biguanide (Polyhexanide)*.

Appearance and smell: clear, colourless and virtually odourless, aqueous solution. **Shelf life:** according to the expiry date, store at room temperature. **Originality:** sterile, originality seal. **Ampoules:** for single use only. **Bottles:** for single patient use. **Manufacturing site:** Made in Germany.

*if used separately, can be considered to be a medicinal product.

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B. Braun Australia Pty Ltd
**B. Braun Medical AG**