



Veterinary Infusion Pump

Proven Excellence for **Premium Quality Medicine**



Brand Recognition & Proven Quality

Our brand reliability and superior quality will win the trust of pet owners and strengthen the reputation of your practice.


Safe to Use

Tube Fatigue Compensation

AniFM I1 can ensure the accuracy of infusion after prolonged or repeated use of the infusion tube

Occlusion Management

AniFM I1 is able to deal with the occlusion caused by tube bending or animal pressure, and automatically resume the infusion after the occlusion resolves, while eliminating post-resolution (unintended) bolus delivery.





Air Bubble Detection

The Air Bubble Detection has 5 detection levels with different sensitivity to detect different sized air bubbles and reduce false alarms

Easy to Operate







Desktop



Pole Clamp

Easy operation with Quick guide on the right panel

One touch flow rate change, no need to Start/Stop repeatedly

AC socket on the left side avoids chewing of the power cord













www.mindrayanimal.com service@mindrayanimal.com

P/N:ENG-AniFM I1-ORA-210285X2P-20241016

